

Sophisticated Data Analysis

FCS Express is packed with tools specifically designed for today's complex flow cytometry data analysis.



FCS Express lets you get your results quickly and easily, without you wasting time with tedious analysis details.

Advanced Gating

Hierarchical and Boolean gating provides complete flexibility in your gating strategy. Use gate linking to analyze in batches or break the link to make adjustments on a file by file basis.

Create your own Statistics

Drag and drop plots onto text boxes to create custom stats tables. Create formulas to define your own stats that automatically update as you change your analysis.



Software Compensation

Adjust the compensation yourself or let *FCS Express* calculate the proper compensation for you.

Advanced Quadrants

Choose between standard or non-rectangular quadrants. Quadrants can automatically adjust based on your data. Link quadrants on different plots or use them as gates.

Powerful Tools for Analysis

Concentration Calculator – calculate the absolute concentration of any population in your sample.

High Throughput Analysis – analyze well based data with heat maps and well gates.

Parameter Math – create new parameters that automatically update during analysis.

Channel Calibration – calibrate your axes and statistics in terms of molecules of fluorochrome.

Histogram Subtraction – automatically calculate the percent positive of your sample.

For features, pricing and your FREE 30-day evaluation, visit:
www.denovosoftware.com



* MS Office is a registered trademark of Microsoft Corporation.